

Liko Stretchers



- OctoStretch is one of Liko's most common stretchers and gets its name from its eight arms
- With Liko's Stretch Leveller, individual adjustment of the centre of gravity for each patient is easy
- The centre of gravity is shifted in the desired direction by simply turning the knobs



Compliance

AS/NZS ISO 10535:2011

Warranty Information

- Frame: 10 years



| Code | Display Name | Safe Working Load | Unit Weight |
|-----------|--|-------------------|-------------|
| LSS395220 | Stretch Leveller - Fits Mod 600 IC and OctoStretch | 200kg | - |
| LSS395230 | Rigor Stretcher | 150kg | 10.3kg |