

Body Assist 42S Knee Support



- An open patella with medial & lateral medical side stays for extra overall support
- Open Patella minimises pressure over the knee cap
- Quality lightweight elastic orthopaedic soft tissue support with additional side supports

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Knee Circumference
FBS284260	BodyAssist - Knee Brace - Elastic Slip-on Support - Sml	305 - 355 mm
FBS284295	BodyAssist - Knee Brace - Elastic Slip-on Support - XXL	460 mm+
FBS284290	BodyAssist - Knee Brace - Elastic Slip-on Support - XL	405 - 460 mm
FBS284270	BodyAssist - Knee Brace - Elastic Slip-on Support - Med	330 - 380 mm
FBS284280	BodyAssist - Knee Brace - Elastic Slip-on Support - Lge	335 - 405 mm