

## Dorset Chair



- A classic and sophisticated chair, the Dorset is perfectly suited to mid-century and modern interiors
- The timeless design features simple back upholstered in suede-look fabric
- Elevated natural timber legs offer a vintage touch
- The padded seat and high back deliver maximum comfort and support
- A refined and versatile chair suited to a lounge or waiting area
- Colour/material to be advised at time of order



### Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Total Height
FUL052130	Dorset Chair	650mm	595mm	870mm