

## TabTimer The Handi Pill Organiser Timer

### TabTimer The Handi Pill Organiser Timer



- Up-to 4 alarms per day
- 28 compartments for storing different medications for the whole week
- Easy to set and use
- When due, alarms sounds for 1 minute
- Then a missed pill reminder, one beep every 20 seconds
- Press any button to stop
- 12 / 24 hour clock
- LCD display with backlight
- Volume level approx 45-50db
- 2 x AAA size battery included
- LCD Viewing Area: 16 x 22mm
- Timer weight 29g without batteries, complete compartments & case 96g
- Timer Dimensions: 28mm x 20mm x 90mm (without 7 days compartments)
- Complete Compartments & Case Dimensions: 186mm x 105mm x 22mm
- Daily compartments: 87mm x 15mm x 18mm
- Individual dose compartments: 20mm x 15mm x 18mm, room for up to 5 panadols
- Requires reasonable dexterity to open compartments

**Code**

**Display Name**

DLG275060 TabTimer The Handi Pill Organiser - 4 Alarm 28 Compartment Pill Box Timer