


## Colonial Weigh Chair Digital



- Lift away armrests and non-slip foot support
- Push handles and 4" castors for portability
- Safety stops on rear castors for parking without rolling
- Changeable plastic seat for easy cleaning
- Internal rechargeable battery
- Seat dimensions: **530 mm width x 350 mm depth**

### Warranty Information

- 1 Year

 Contract approved

### Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
LSS377704	Colonial Weigh Chair	840 mm	665 mm	940 mm	250 kg