

Leg Lifter



- Enables the user to move a stiff or immobile leg
- Upper loop fits around the hand and wrist and the lower stiffened loop slides over the foot to give control when lifting the leg
- Stiffened stem assists with manoeuvring leg
- Ideal assistance for getting into or out of a bed, wheelchair or using a footstool etc.
- Weight: 80 g

Warranty Information

• 1 Year

0			
\odot	Contracted	DVA	produc

⊘Contract approved

ndis

Item Code:

05_090900046_0103_1_2 Click **here** for more information

Compliance

No applicable compliance standard

Code	Display Name	Total Length
DI D269390	Leg Lifter	864 mm