

## Leg Lifter



- Enables the user to move a stiff or immobile leg
- Upper loop fits around the hand and wrist and the lower stiffened loop slides over the foot to give control when lifting the leg
- Stiffened stem assists with manoeuvring leg
- Ideal assistance for getting into or out of a bed, wheelchair or using a footstool etc.
- Weight: 80 g

### Warranty Information

- 1 Year

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Dressing and undressing

**Category:** Low risk

Click [here](#) for more information

### NDIS

**Item Code:**

05\_090900046\_0103\_1\_2

Click [here](#) for more information

### Compliance

No applicable compliance standard

| Code      | Display Name | Total Length | Unit Weight |
|-----------|--------------|--------------|-------------|
| DLD269390 | Leg Lifter   | 864 mm       | 80 g        |