

Leg Lifter



- Enables the user to move a stiff or immobile leg
- Upper loop fits around the hand and wrist and the lower stiffened loop slides over the foot to give control when lifting the leg
- Stiffened stem assists with manoeuvring leg
- Ideal assistance for getting into or out of a bed, wheelchair or using a footstool etc.
- Weight: 80 g

Warranty Information

- 1 Year

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Dressing and undressing

Category: Low risk

Click [here](#) for more information

NDIS

Item Code:

05_090900046_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Total Length
DLD269390	Leg Lifter	864 mm