

## Leg Lifter



- Enables the user to move a stiff or immobile leg
- Upper loop fits around the hand and wrist and the lower stiffened loop slides over the foot to give control when lifting the leg
- Stiffened stem assists with manoeuvring leg
- Ideal assistance for getting into or out of a bed, wheelchair or using a footstool etc.
- Weight: 80 g

### Warranty Information

- 1 Year

- ✓ Contracted DVA product
- ✓ Contract approved

### ndis

#### Item Code:

05\_090900046\_0103\_1\_2

Click [here](#) for more information

### Compliance

No applicable compliance standard

Code	Display Name	RRP	Box Quantity	Total Length
DLD269390	Leg Lifter	39.2	1	864 mm