

Dane Chair



- The distinct back of this Dane Chair perfectly captures a retro aesthetic
- High back and padded cushion delivers maximum comfort and support
- Retro design is complemented by a mid-century style timber base
- Upholstered in a soft and durable fabric with a subtle, on-trend print
- A unique occasional chair for a living or waiting area
- Fabric colour/material to be advised at time of order



Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Total Height
FUL052120	Dane Chair	585mm	710mm	920mm