

Oppo 4039 - Motion Control ROM - Open Knee Splint



- Provides stability to an injured or unstable knee by allowing normal knee joint motion with proper settings or stops within a selected range for use following knee surgery, early cast removal or post-operative rehabilitation
- Settings from 0 to 120 degrees flexion and between 0 to 75 degrees extension in 7.5 degree increments
- Hook and loop closures for uniform compression and support
- Application: post surgical use, collateral ligament injury, cruciate ligament injury, knee instability
- Available in three lengths: 46cm, 56cm and 66cm

 Contract approved

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Total Length	Colour	Knee Circumference
FBS301536	Oppo 4039 - Motion Control ROM - Open Knee Splint - 460 mm	460 mm	Blue	432 - 711 mm
FBS30153622	Oppo 4039 - Motion Control ROM - Open Knee Splint - 560 mm	560 mm	Blue	432 - 711 mm
FBS30153626	Oppo 4039 - Motion Control ROM - Open Knee Splint - 660 mm	660 mm	Blue	432 - 711 mm