

Oppo 1031 Knee Stabiliser Hinged



- Neoprene circle pad around the patella provides padded protection and stabilization to prevent patella displacement.
- Two medial and lateral stays provide additional support and stability for the knee.
- Two hook and loop straps insure comfortable fit and stabilization during activities.
- Helps relieve pain from tendinitis and arthritis.
- Retains body heat to warm knee joint and increases blood circulation.
- For accurate sizing measure around mid section of the knee, while standing.
- Application:
 - Patellar subluxation/tilting
 - Collateral ligament injury
 - Cruciate ligament sprain
 - Moderate sprain/ strain

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Knee Circumference
FBS301528	Oppo 1031 Knee Stabiliser Hinged - S	333 - 375 mm
FBS301528L	Oppo 1031 Knee Stabiliser Hinged - L	406 - 451 mm
FBS301528M	Oppo 1031 Knee Stabiliser Hinged - M	375 - 406 mm
FBS301528XL	Oppo 1031 Knee Stabiliser Hinged - XL	451 - 502 mm