

Jobst Bella Strong - Gauntlet (20 - 30 mmHg)



- Ideal for managing moderate upper extremity lymphoedema
- Donning device included with every sleeve for easy donning
- Open knit structure allows for improved air permeability which keeps skin dry
- Two-way stretch makes the arm sleeve comfortable without compromising medical efficacy
- Circular knit construction
- Cotton-like feel; feels soft to the skin
- Dotted, 2-inch silicone band keeps sleeve firmly in place

- ✓ Contracted DVA product
- ✓ Contract approved



Code	Display Name	Wrist Circumference
PTS642534	Jobst Bella Strong - Ready To Wear - Gauntlet 20-30mmHg - Sz 1	140 - 180 mm
PTS642535	Jobst Bella Strong - Ready To Wear - Gauntlet 20-30mmHg - Sz 2	140 - 180 mm
PTS642536	Jobst Bella Strong - Ready To Wear - Gauntlet 20-30mmHg - Sz 3	180 - 210 mm
PTS642537	Jobst Bella Strong - Ready To Wear - Gauntlet 20-30mmHg - Sz 4	140 - 180 mm
PTS642538	Jobst Bella Strong - Ready To Wear - Gauntlet 20-30mmHg - Sz 5	180 - 210 mm
PTS642539	Jobst Bella Strong - Ready To Wear - Gauntlet 20-30mmHg - Sz 6	160 - 190 mm
PTS642540	Jobst Bella Strong - Ready To Wear - Gauntlet 20-30mmHg - Sz 7	190 - 220 mm

Code	Display Name	Wrist Circumference
------	--------------	---------------------

PTS642541 Jobst Bella Strong - Ready To Wear - Gauntlet 20-30mmHg - Sz 8 190 - 220 mm