

## Aspire Female Urinal Cygnet



Aspire Female Urinal Cygnet are fully autoclavable and made from premium polypropylene. They will last for many years in any healthcare environment subject to the rigours of cleaning, sanitation and infection control requirements.

- Heavy duty and reusable
- Designed for use in healthcare environments such as surgeries, acute care, general ward, medical centres, and aged care facilities
- Autoclavable polypropylene that can be cleaned using autoclaves and sanitisers
- Pastel green colour for easy identification

### Warranty Information

- 1 Year



- ✔ Contracted DVA product
- ✔ Contract approved

### Support @ Home

**Function:** Collecting urine and faeces

**Category:** Prescribed

Click [here](#) for more information

### NDIS

**Item Code:**

05\_098800044\_0103\_1\_2

**Toilet Aids – Composite:**

05\_414\_0103\_1\_2

Click [here](#) for more information

### Compliance

AS/NZS 3973:2009

Code	Display Name
BTT148800	Female Urinal Cygnet

