

## Aspire ComfiMotion Fixed Companion Bed Package



The Aspire ComfiMotion Fixed Companion Bed Package is the perfect blend of style and practicality, designed to complement the ComfiMotion Adjustable Bed range. This package includes the ComfiMotion Fixed Companion Bed together with the Aspire ComfiMotion Allure or Chesterfield Free Standing Headboard as standard, creating a modern, coordinated look that enhances both comfort and bedroom design.

 Contract approved

- Matching Long Single bed without any adjustability designed to sit flush when positioned together with the Aspire ComfiMotion Care Bed and Aspire ComfiMotion Activ Care Bed.
- It facilitates a side-by-side sleeping arrangement at a fixed height of 370 mm, without compromising on the users care and assistance requirements.



### Warranty Information

- 10 Years - Frame
- 1 Year - Fabric

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Colour
BEB046800V2	Aspire ComfiMotion Fixed Companion Bed Package - Chesterfield - Long Single	2030 mm	920 mm	370 mm	180 kg	Portobello Grey
BEB046800V1	Aspire ComfiMotion Fixed Companion Bed Package - Allure - Long Single	2030 mm	920 mm	370 mm	180 kg	Portobello Grey