

Aspire Adjustable Day Chair



Jaspire

- A high back orthopaedic chair that offers adjustability, comfort and support
- Depth Adjustable Ensures changing user needs are accommodated
- Contoured Backrest High back design provides increased postural support
- Padded Seat and Armrests Generous support surfaces for better pressure distribution and sitting comfort
- Padded Headrest Cushion Provides additional comfort for long periods of sitting
- Sturdy Design High safe working load of 160 kg, yet lightweight for easy manoeuvrability
- Transport Wheels Allows for the chair to be moved with ease by carers
- Height Adjustable Easy push button adjustment for comfortable positioning and changing needs
- Non-Slip Tips Provide stability whilst sitting or transferring
- Powder Coated Steel Frame Easy to clean, strong and durable

Warranty Information

• 1 Year

- Ocontracted DVA product
- Ocontract approved

ndis

Item Code:

05_188800051_0103_1_2 Click **here** for more information

Compliance

AS/NZS 4688.2:2000; AS/NZS 4688.3:2000

Code	Display Name	Seat Depth	Seat Width	Total Height	Seat To Floor Height	Safe Working Load	Unit Weight
CHP208025	Aspire Adjustable Day Chair - Latte Vinyl	475 - 590 mm	510 - 610 mm	1070 - 1220 mm	410 - 530 mm	160 kg	26 kg
CHP208030	Aspire Adjustable Day Chair - Ink Vinyl	475 - 590 mm	510 - 610 mm	1070 - 1220 mm	410 - 530 mm	160 kg	26 kg