

Aspire ComfiMotion Relax Bed Package - Chesterfield



The Aspire ComfiMotion Relax Bed Package offers the perfect balance of comfort, support and contemporary style for the home. Complete with the Chesterfield Free Standing Headboard as standard, this package enhances the Relax Bed's adjustable features with a sleek, modern finish - delivering a bed that looks as good as it feels.

The Relax Bed has a fixed base and combines luxurious finishes with multiple positioning and therapeutic functions to support users and carers alike.

- **Night Assist:** Luminous underbed halo for safe entry and exit throughout the night.
- **Massage:** Relaxing vibration modes for personalised neck, shoulder and back relief.
- **Box-on-Box Design:** Facilitates easier entry and exit, allowing true flush positioning of beds in companion or dual adjustable bed configurations.
- **German Engineered Motors:** For smooth, quiet adjustment.
- **USB Charge:** Connect and charge your devices conveniently.
- **Head Tilt:** Relaxes neck muscles, providing optimal support whilst reading or watching TV.
- **Back Raise:** Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea.
- **Leg Raise:** Promotes circulation, relieving tired, swollen and/or painful limbs.
- **Zero Gravity:** Helps to relax the circulatory system and muscles of the lower back.
- **Memory Function:** Allows the user to set and return to their preferred positions at any time with one click of a button.
- Total height adjustable in increments of 320, 370, 380 & 430 mm



Warranty Information

- 10 Years Frame
- 2 Years Electronics
- 1 Year Handset, Cables & Fabric
- 1 Year Headboard

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour
BEB046910V2	Aspire ComfiMotion Relax Bed Package - Chesterfield - King Single (inc. Freestanding Headboard)	2030 mm	1070 mm	320 - 430 mm	180 kg	58 kg	Portobello Grey
BEB046920V2	Aspire ComfiMotion Relax Bed Package - Chesterfield - Long Double (inc. Freestanding Headboard)	2030 mm	1370 mm	320 - 430 mm	360 kg	62 kg	Portobello Grey
BEB046930V2	Aspire ComfiMotion Relax Bed Package - Chesterfield - Queen (inc. Freestanding Headboard)	2030 mm	1520 mm	320 - 430 mm	360 kg	70 kg	Portobello Grey
BEB046900V2	Aspire ComfiMotion Relax Bed Package - Chesterfield - Long Single (inc. Freestanding Headboard)	2030 mm	920 mm	320 - 430 mm	180 kg	55 kg	Portobello Grey