

Oppo 1132 Contour Knee Support



- Unique contour shape to snugly fit the knee joint to prevent slipping due to excessive movement.
- The adjustable design provides easy application and maximum support.
- Breathable CoolPreene material lined with coolmax behind the knee for optimal comfort.
- Application:
 - Osteoarthritis
 - Knee pain
 - Knee instability

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Knee Circumference
FBS286040	Oppo 1132 Contour Knee Support - XXL	505 - 545 mm
FBS286045	Oppo 1132 Contour Knee Support - XXXL	545 - 585 mm
FBS301505	Oppo 1132 Contour Knee Support - S	330 - 395 mm
FBS301505L	Oppo 1132 Contour Knee Support - L	430 - 470 mm
FBS301505M	Oppo 1132 Contour Knee Support - M	395 - 430 mm
FBS301505XL	Oppo 1132 Contour Knee Support - XL	470 - 505 mm