

## Cura Alarms - ActiveCare Monitor



**cura**<sup>1</sup>

The ActiveCare Monitor is a sophisticated wireless nurse call alarm monitor that can be setup quickly and easily. It supports pads/mats, call, door type alarms, and signals to pager with caller ID and device type displayed.



- Hassle-free installation, no wiring required.
- Wireless caller display alarm (60 devices).
- Integral signal repeater transmits alarm signals to pager.
- Mains powered with battery backup (plug pack and 3 x AA batteries are included).
- Mains power failure or interruption notified to paired pagers.
- Multiple call display - review calls from the last 15 mins.
- Separate indication for DoorAlert.
- Selectable 'signal lost' tracking of bed, chair, and mat devices.
- Device type indication during active alarms.

### Warranty Information

- 1 Year

Code	Display Name
BEA007774	Cura - ActiveCare Monitor