

Freedom Oval Tube Shower Stool - HBA408



freedom
HEALTHCARE

- Moulded seat, perforated for drainage
- Height adjustable legs. Adjustment pins are spaced at 25 mm increments on three extension legs and 12.5 mm increments on the fourth to maximise stability
- S-bend shaped frame, provides maximum space between the armrests
- Powder coated aluminium frame
- Slip resistant suction rubber tips
- Purchased separately: Extension Legs (GSP650100), Armrest Pad (GSP650460)

✓ Contracted DVA product

✓ Contract approved

Compliance

AS/NZS 3973:2009

Warranty Information

- 1 Year

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height
BTS150140	Freedom Oval Tube Shower Stool - 200 kg	242	1	520 mm	610 mm	670 - 820 mm	200 kg	4.2 kg	290 mm	400 mm	480 - 630 mm