

Freedom Oval Tube Shower Stool - HBA408



- Moulded seat, perforated for drainage
- Height adjustable legs. Adjustment pins are spaced at 25 mm increments on three extension legs and 12.5 mm increments on the fourth to maximise stability
- S-bend shaped frame, provides maximum space between the armrests
- Powder coated aluminium frame
- Slip resistant suction rubber tips
- Purchased separately: Extension Legs (GSP650100), Armrest Pad (GSP650460)

Warranty Information

- 1 Year

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height
BTS150140	Freedom Oval Tube Shower Stool - 200 kg	520 mm	610 mm	670 - 820 mm	200 kg	4.2 kg	290 mm	400 mm	480 - 630 mm