

DermaSaver Knuckle Protector



- Designed for comfort and effectiveness to protect fragile skin on the knuckles
- Constructed of two layers of soft, lightweight MicroSpring textile, which allows air-flow and moisture transfer
- Soft and loose with ample air-flow reaching the skin
- Anti-microbial - will not support growth of bacteria, latex free
- Available in 4 sizes to suit small to extra large hands
 - Medium - Knuckle Circumference < 20cm
 - Large - Knuckle Circumference > 20cm

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Code	Display Name
PTA533301	DermaSaver - Knuckle Protector - Sml
PTA533301L	DermaSaver - Knuckle Protector - Lge