

Aspire ComfiMotion Care Bed



- Night Assist: Luminous underbed halo for safe entry and exit throughout the night
- Press and hold safety handset: Intuitive, large button interface enhances usability
- Enclosed frame design: Ensures optimal positioning when getting into the adjustable bed and safely covers the moving base platform
- Safe-Lowering Technology: Ensures that back, leg and Hi-Lo functions only lower under the force of gravity reducing potential for accidents or injuries
- Box-on-box design facilitates easier entry and exit and allows true flush positioning of beds in companion or dual adjustable bed configurations.
- German engineered motors for smooth, quiet adjustment
- *Total Height with:
 - Low Profile Fixed Legs 280 – 580 mm (available upon request)
 - 75mm Standard Castors 350 - 650 mm
 - 125mm Accessibility Castors 400 - 700 mm (optional)Mattress choice may increase or decrease overall bed height
- *Mattress not included

Prescribing Therapists: Click here to download our
[Clinical Assistance Guide](#)

Warranty Information

- 10 Years Limited Warranty

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click [here](#) for more information

Compliance

IEC 60601-2-52:2009

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
BEB046720	Aspire ComfiMotion Care Bed - Portobello Grey - Long Single (Headboard Sep)	2030 mm	920 mm	280 - 580 mm	180 kg
BEB046730	Aspire ComfiMotion Care Bed - Portobello Grey - King Single (Headboard Sep)	2030 mm	1070 mm	280 - 580 mm	180 kg