

BodyAssist Elbow Sleeve



- Provides a basic level of compression and orthopaedic support

✔ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Colour	Elbow Circumference
FBS284640	BodyAssist 95 Elbow Sleeve - Small	Beige	200 - 250 mm
FBS284650	BodyAssist 95 Elbow Sleeve - Medium	Beige	250 - 280 mm
FBS284660	BodyAssist 95 Elbow Sleeve - Large	Beige	280 - 300 mm
FBS284670	BodyAssist 95 Elbow Sleeve - XLarge	Beige	300+ mm