


K Care Pressure Reducing Commode Seat



- Special foam design, engineered to improve the spread of pressure across the seat
- Compliments the user's pressure management routine
- Waterproof but breathable antimicrobial four-way stretch Dartex material to optimises user comfort
- Independently tested and proven to be effective in reducing pressure levels
- Seat height: 65 mm

Warranty Information

- 1 Year



Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

Code	Display Name	Seat Depth	Seat Width
BTC082100	Closed Front Commode Seat - Pressure Reducing - 440 mm Wide	430 mm	440 mm
BTC082000	Open Front Commode Seat - Pressure Reducing - 440 mm Wide	430 mm	440 mm