


## TabTimer Talking Electronic Body Weight Scale

### TabTimer Talking Electronic Body Weight Scale



- Talking audible and digital weight results
- Equipped with a high precise strain gauge sensor system
- Tempered safety glass platform
- Auto zero resetting and auto power off
- Low power and over-load indicator
- Unit button to convert into kg/lb/st.
- Maximum weight 176kg
- Dimensions: 280 x 280 x 20mm
- 2 x AAA Batteries



Code	Display Name
DLG280000	TabTimer Talking Electronic Body Weight Scale