

Medidisc Transfer Turn Table



- A floor turntable for use when transferring people who can stand but have difficulty turning
- The user is pivoted around with the assistance of an attendant
- The MediDisc has slip resistant studs on the base, a flexible rim and a slot handle for transportation
- Available in one standard size only - 380 mm Diameter

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Changing body position

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Safe Working Load	Unit Weight
TSE672940	Medidisc Transfer Turn Table	115 kg	1 kg