

Oppo 1036 Knee Support Multi Orthosis



- Functional 3/16' breathable neoprene knee brace.
- Polycentric hinge with flexion and extension pin-stops for controlled range of motion.
- Lightweight, open patella design.
- Increases stability and provides support and protection.
- Provides greatest comfort even in the hollow of the knee.
- For patients with restricted grip force and mobility.
- For accurate sizing measure around mid section of the knee, while standing.
- Application:
 - Post surgical use
 - Patella subluxation / tilting
 - Collateral ligament injury
 - Cruciate ligament sprain
 - Moderate sprain / strain

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Knee Circumference
FBS286030	Oppo 1036 Knee Support Multi Orthosis - XXL	502 - 546 mm
FBS301534	Oppo 1036 Knee Support Multi Orthosis - S	333 - 375 mm
FBS301534L	Oppo 1036 Knee Support Multi Orthosis - L	406 - 451 mm
FBS301534M	Oppo 1036 Knee Support Multi Orthosis - M	375 - 406 mm
FBS301534XL	Oppo 1036 Knee Support Multi Orthosis - XL	451 - 502 mm