

ETAC Cloo Fixed Height Toilet Seat Raiser with Arms



- The Cloo from Etac is a height adjustable toilet seat raiser with tiltable seat and detachable armrests for added support
- The arm supports fold up and away, or can be fixed in a folded down position for extra security if needed
- Both armrests are detachable, Cloo can also be used with only one arm support if space around the toilet is limited
- The seat can be easily adjusted to three different heights: low (60 mm), medium (100 mm) or high (140 mm) for maximum versatility
- Simple, elegant toilet seat raiser
- Detachable armrests either side of seat
- Mount to your existing toilet three height settings
- Easy to adjust height while mounted
- Can be set with forward facing tilt
- Detachable seat for easy cleaning
- Tool-less mounting to any toilet bowl

 Contract approved

Support @ Home

Function: Toileting

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Total Width	Total Height	Safe Working Load	Unit Weight	Seat Width
BTT162030	ETAC Cloo Fixed Height Toilet Seat Raiser with Arms	600 mm	60 - 140 mm	150 kg	4.8 kg	400 mm