

DermaSaver Shin-Knee Tube



- Two layers of soft, lightweight MicroSpring textile that protect fragile skin from pressure, friction and minor traumas that result in tears, bruising and abrasions
- Designed especially for residents in wheelchairs to protect both the shin & the knee
- The narrower end of the tube gently rests on the ankle to keep it from falling off

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Calf/Arm Circumference
PTA533323	DermaSaver - Shin-Knee Tube - X Sml	18-23cm
PTA533323S	DermaSaver - Shin-Knee Tube - Sml	23-30cm
PTA533323M	DermaSaver - Shin-Knee Tube - Med	30-38cm
PTA533323L	DermaSaver - Shin-Knee Tube - Lge	38-48cm
PTA533323XL	DermaSaver - Shin-Knee Tube - XL	48-58cm