

Scooter Oxygen Bottle Holder



PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

- Scooter mounted oxygen bottle holder
- Carries a 406L oxygen bottle
- Fits easily into the scooter tow bar on the rear
- Can be used in conjunction with the cane holder mount
- If by itself, requires the "L" Piece Accessory (SCA668310)

Warranty Information

- 1 Year

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Diverse motorised vehicles

Category: Prescribed

Click [here](#) for more information

Compliance

AS3695.2:2013

Code	Display Name
SCA668330	Oxygen Bottle Holder