

Scooter Oxygen Bottle Holder



PRIDE
MOBILITY
LIVE YOUR BEST®

- Scooter mounted oxygen bottle holder
- Carries a 406L oxygen bottle
- Fits easily into the scooter tow bar on the rear
- Can be used in conjunction with the cane holder mount
- If by itself, requires the "L" Piece Accessory (SCA668310)

Warranty Information

- 1 Year

✓ Contracted DVA product

✓ Contract approved

Compliance

AS3695.2:2013

Code	Display Name	RRP	Box Quantity
SCA668330	Oxygen Bottle Holder	186	1