

## Aspire LIFESTYLE Community Bed



Adjustable and stylish Homecare bed designed for easy transportation with a modern look suited to the Homecare Environment.

- Warm timber finishes paired with a modern dark frame, the Lifestyle bed not only looks sophisticated in the home environment, it features five profiling functions that cater to a broad spectrum of comfort and positioning needs for the ultimate in home care provision
- **Quick setup time:** Tool free with no requirement to remove electronics during assembly and set-up, ensuring quick and easy readiness upon delivery or for transportation and storage. For someone working alone it should take a maximum of 8 minutes to assemble a Lifestyle bed!
- **Additional sizes:** available in both single and king single widths to suit most customer needs
- **Homecare upgrade options:** Timber laminate head and foot board shrouds and full-length timber rails (as pictured) can be fitted for an appealing homely look (accessories purchased separately click [here](#))
- **Supplied Transport-Ready:** The Lifestyle beds come pre-assembled on transport brackets, so you can simply roll it out of the carton on its wheels, ready to deliver and set-up!
- **Transportable Design:** Tool free with no requirement to remove electronics during assembly and set-up, ensuring quick and easy readiness upon delivery or for transportation and storage
- **Leading Quality Frame & Electronics:** Strong durable steel frame with Linak electronics and handset allows for market leading warranty provision
- **Four Total-Lock Castors:** Safely locate the bed and prevent castor roll and directional rotation
- **Highly Adjustable** Four electric profiling functions to suit all comfort and positioning needs

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Beds and bed equipment

**Category:** Prescribed

Click [here](#) for more information

### Accessories

Click on the cart icon to display the range of optional accessories

- The robust and easily assembled frame makes for the perfect rental fleet bed, built to withstand the rigors of regular transportation and assembly.
- Full Safe Working Load: 220 kg, Maximum Patient Weight 180 kg

## Warranty Information

- Frame: 3 Years
- Electronics: 2 Years
- Handset: 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Colour	Number of Functions
BEB046030	Aspire LIFESTYLE Community Bed w/Transport Bracket - Single	2220 - 2370 mm	995 mm	245 - 645 mm	220 kg	Gunmetal Grey	5
BEB046040	Aspire LIFESTYLE Community Bed w/Transport Bracket - King Single	2220 - 2370 mm	1155 mm	245 - 645 mm	220 kg	Gunmetal Grey	5