

## Aspire Slipper Bed Pan



Aspire Bed Pan Slipper's are fully autoclavable and made from premium polypropylene. They will last for many years in any healthcare environment subject to the rigours of cleaning, sanitation and infection control requirements.

- Heavy duty and reusable
- Designed for use in healthcare environments such as surgeries, acute care, general ward, medical centres, and aged care facilities
- Autoclavable polypropylene that can be cleaned using autoclaves and sanitisers
- Pastel green colour for easy identification

### Warranty Information

- 1 Year



- ✓ Contracted DVA product
- ✓ Contract approved

### ndis

#### Item Code:

05\_098800044\_0103\_1\_2

#### Toilet Aids – Composite:

05\_414\_0103\_1\_2

Click [here](#) for more information

### Compliance

AS/NZS 3973:2009

Code	Display Name	RRP	Box Quantity
BTT148000	Bed Pan Slipper - Small	43.4	1
BTT161500	Bed Pan Slipper - Large	47.4	1