

Seahorse Plus Hygiene Chair



- Gas assisted 40° Tilt-In-Space to enable showering, hair washing and toileting
- Tilt locking safety mechanism - simple push button locking tilting mechanism for ease of operation and security
- 100 mm height adjustment
- Height adjustable flip-up footrests with foot straps
- Height and depth adjustable headrest
- Removable armrests, flip-up anterior supports, pommel and potty
- Contoured seat with soft padding provides additional comfort and aids relaxation
- Designed for use over most standard toilets
- The seat can be raised or lowered for optimum use
- Large castors with brakes for ease of movement and maximum security
- Removable two-piece liners: Small/Medium (2 sizes); Large (3 sizes)
- Contoured headrest
- Shoulder harness
- Leg Strap



| Code | Display Name | Total Width | Safe Working Load | Unit Weight | Seat Depth | Seat Width | Seat To Floor Height | Seat To Footrest |
|-----------|------------------------|-------------|-------------------|-------------|------------|--------------|----------------------|------------------|
| PED479650 | Seahorse Plus - Small | 532 mm | 63.5 kg | 18.7 kg | 374 mm | 332 - 365 mm | 521 - 629 mm | 521 - 629 mm |
| PED479660 | Seahorse Plus - Medium | 586 mm | 88.9 kg | 22.3 kg | 470 mm | 378 - 423 mm | 501 - 609 mm | 501 - 609 mm |
| PED479665 | Seahorse Plus - Large | 692 mm | 102 kg | 24.8 kg | 488 mm | 477 - 530 mm | 501 - 609 mm | 501 - 609 mm |