

Seahorse Plus Hygiene Chair



- Gas assisted 40° Tilt-In-Space to enable showering, hair washing and toileting
- Tilt locking safety mechanism - simple push button locking tilting mechanism for ease of operation and security
- 100 mm height adjustment
- Height adjustable flip-up footrests with foot straps
- Height and depth adjustable headrest
- Removable armrests, flip-up anterior supports, pommel and potty
- Contoured seat with soft padding provides additional comfort and aids relaxation
- Designed for use over most standard toilets
- The seat can be raised or lowered for optimum use
- Large castors with brakes for ease of movement and maximum security
- Removable two-piece liners: Small/Medium (2 sizes); Large (3 sizes)
- Contoured headrest
- Shoulder harness
- Leg Strap

✓ Contract approved

Support @ Home

Function: Toileting

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Total Width	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
PED479650	Seahorse Plus - Small	532 mm	63.5 kg	18.7 kg	374 mm	332 - 365 mm	521 - 629 mm	521 - 629 mm
PED479660	Seahorse Plus - Medium	586 mm	88.9 kg	22.3 kg	470 mm	378 - 423 mm	501 - 609 mm	501 - 609 mm
PED479665	Seahorse Plus - Large	692 mm	102 kg	24.8 kg	488 mm	477 - 530 mm	501 - 609 mm	501 - 609 mm