

## Everyday Seat

**LECKEY®**

### Paediatric Therapy Chair

#### Optimal Posture + Stability = Everyday

- Why - Additional pelvic support benefits those with moderate postural control needs because the pelvis provides the key to overall postural stability
- How - The Everyday Activity Seat offers exceptional pelvic positioning with a four-point pelvic harness which is easily adjusted from the centre buckle. This is complemented by the seat's simple range of trunk, head, leg and foot positioning adjustments. In keeping with the needs of this user group, the Everyday Activity Seat provides optimal postural stability

#### Optimal Comfort = Everyday

- Why - It is not possible to sustain an optimal posture without comfort because discomfort disturbs posture
- How - The Everyday Activity Seat is designed with the comfort of the young person in mind in every aspect. From the soft, breathable covers, to the moulded armrests, and the seat base cover which tucks underneath the seat base, the contoured back and the ramped seat cushion, the Everyday Activity Seat ensures the comfort (and therefore postural stability) of the client

#### Optimal Function = Everyday

- Why - Function (the ability to do everyday things) depends on the combination of postural stability and comfort
- How - With optimal posture, stability and comfort achieved, the client will be able to use their arms and hands more effectively and independently. Along

**ndis****Item Code:**

05\_180939187\_0103\_1\_2

Click [here](#) for more information

with the Everyday Activity Seat's height adjustability from floor to table (hi-low chassis), activities such as playing, eating, reading, writing or chatting with friends become easier and less tiring

## Warranty Information

- 3 Years

Code	Display Name	Seat Depth	Seat Width	Backrest Height	Handle Height	Safe Working Load	Colour
PED490015	Everyday Seat - Size 2	230 - 380 mm	260 - 380 mm	400 - 520 mm	200 - 300 mm	50 kg	-
PED490020	Everyday Seat - Black - Size 2	230 - 380 mm	260 - 380 mm	400 - 520 mm	120 - 150 mm	50 kg	Black