

Glove Applicator



- Accessory for donning and doffing medical compression garments
- Specially designed to provide an optimal grip on the fabric and help medical compression stockings to slide over the foot and up the leg
- Preserves the original medical compression by avoiding pulling hard at the top of your stockings and ensuring a longer life of the garment
- Protects the fabric of stockings from being damaged by nails or jewellery
- Contains Latex

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Dressing and undressing

Category: Low risk

Click [here](#) for more information

Code	Display Name
PTS643440	Stocking Applicator - Gloves - Small
PTS643430	Stocking Applicator - Gloves - Medium
PTS643420	Stocking Applicator - Gloves - Large
PTS643450	Stocking Applicator - Gloves - Xlarge