

Ottobock Evolight



- Lightweight wheelchair cushion.
- Cube structured foam core with horizontally integrated air channels for improved air circulation.
- Two layers of spacer fabric enable moisture wicking.
- The cushion attaches to a wheelchair with a hook and loop strap.
- Cover is removable. Suitable for people at low risk of developing a pressure injury.

- ✓ Contracted DVA product
- ✓ Contract approved

ottobock.

Code	Display Name
PTC551000	Evolight - Otto Bock - Cover (Specify Size)
PTC550275	Evolight Cushion - 16x16in
PTC550300	Evolight Cushion - 22x19in
PTC550310	Evolight Cushion - 23x18in