

## Ottobock Evolight



- Lightweight wheelchair cushion.
- Cube structured foam core with horizontally integrated air channels for improved air circulation.
- Two layers of spacer fabric enable moisture wicking.
- The cushion attaches to a wheelchair with a hook and loop strap.
- Cover is removable. Suitable for people at low risk of developing a pressure injury.

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Manage tissue integrity

**Category:** Prescribed

Click [here](#) for more information

**ottobock.**

Code	Display Name	Total Length	Total Width
PTC551000	Evolight - Otto Bock - Cover (Specify Size)	-	-
PTC550275	Evolight Cushion - 16x16in	16 in	16 in
PTC550300	Evolight Cushion - 22x19in	19 in	22 in
PTC550310	Evolight Cushion - 23x18in	18 in	23 in