

MEDI R.141 RHIZOMED SOFT



- Effective stabilization of the metacarpal phalangeal and carpometacarpal joint of the thumb
- Integrated stabilization provides reliable immobilization
- Two adjustable thumb straps to fit all handshapes, providing increased stability
- The splint's semi-open shape adapts to various degrees of swelling of the thumb
- The thumb-index pinch grip and prehension grasp are maintained
- Large opening for ease of application and removal, with as little pain as possible
- Easy to wear, breathable fabric
- Lightweight
- Available in black Application:
- Carpometacarpal arthritis
- Ligament injury (e.g. skier's thumb)
- Inflammation of the thumb joint

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name
FBS307400	MEDI R.141 RHIZOMED SOFT - SIZE I
FBS307410	MEDI R.141 RHIZOMED SOFT - SIZE II
FBS307420	MEDI R.200 WRIST SUPPORT - LEFT - SML
FBS307430	MEDI R.200 WRIST SUPPORT - LEFT - MED

Code	Display Name
FBS307440	MEDI R.200 WRIST SUPPORT - LEFT - LGE
FBS307450	MEDI R.200 WRIST SUPPORT - LEFT - XL
FBS307460	MEDI R.200 WRIST SUPPORT - RIGHT - SML
FBS307470	MEDI R.200 WRIST SUPPORT - RIGHT - MED
FBS307480	MEDI R.200 WRIST SUPPORT - RIGHT - LGE
FBS307490	MEDI R.200 WRIST SUPPORT - RIGHT - XL