

Aspire Bath Transfer Bench - Bariatric



- Reinforced design for the bariatric user with a maximum weight of 227 kg
- It has a textured moulded plastic seat and backrest
- Fully height adjustable with a side arm for secure transfers
- Adjustable backrest to suit right or left facing tubs
- Suction cup foot pieces secure the bench to the tub floor
- Lightweight aluminium frame construction
- **Aspire Extension Legs - [Set of 4](#)**

Each set includes:

- 2 x suction tips (for use on smooth surfaces)
- 2 x standard tips (for general use)

✓ Contracted DVA product

✓ Contract approved

ndis

Item Code:

05_123121160_0103_1_2

Shower Support –

Composite:05_415_0103_1_2

Click [here](#) for more information

Warranty Information

- 1 Year



Compliance

AS/NZS 3973:2009

Code	Display Name	RRP	Box Quantity	Seat Depth	Seat Width	Seat To Floor Height	Safe Working Load	Unit Weight
BTS105300	Aspire Bath Transfer Bench - Bariatric	338	1	405 mm	670 mm	485-595 mm	227 kg	7.2 kg
BTS108312	Extension Legs - 100mm - w/2 Suction & 2 Std Tips - Suits BTS105300	60	1	-	-	-	-	-