

## Aspire Bath Transfer Bench - Bariatric



- Reinforced design for the bariatric user with a maximum weight of 227 kg
- It has a textured moulded plastic seat and backrest
- Fully height adjustable with a side arm for secure transfers
- Adjustable backrest to suit right or left facing tubs
- Suction cup foot pieces secure the bench to the tub floor
- Lightweight aluminium frame construction
- **Aspire Extension Legs - [Set of 4](#)**  
Each set includes:
  - 2 x suction tips (for use on smooth surfaces)
  - 2 x standard tips (for general use)

- ✓ Contracted DVA product
- ✓ Contract approved

### Support @ Home

**Function:** Washing, bathing and showering

**Category:** Prescribed

Click [here](#) for more information

### Warranty Information

- 1 Year

### NDIS

**Item Code:**

05\_123121160\_0103\_1\_2

**Shower Support –**

**Composite:**05\_415\_0103\_1\_2

Click [here](#) for more information

### Compliance

AS/NZS 3973:2009



Code	Display Name	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height
BTS105300	Aspire Bath Transfer Bench - Bariatric	227 kg	7.2 kg	405 mm	670 mm	485-595 mm

Code	Display Name	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height
BTS108312	Extension Legs - 100mm - w/2 Suction & 2 Std Tips - Suits BTS105300	-	-	-	-	-