

## Aspire Bath Transfer Bench - Bariatric





- Reinforced design for the bariatric user with a maximum weight of 227 kg
- It has a textured moulded plastic seat and backrest
- Fully height adjustable with a side arm for secure transfers
- Adjustable backrest to suit right or left facing tubs
- Suction cup foot pieces secure the bench to the tub floor
- Lightweight aluminium frame construction

### Warranty Information

- 1 Year



 Contracted DVA product

 Contract approved

**ndis**

**Item Code:**

05\_123121160\_0103\_1\_2

**Shower Support –**

**Composite:**05\_415\_0103\_1\_2

Click [here](#) for more information

**Compliance**

AS/NZS 3973:2009

Code	Display Name	Seat Depth	Seat Width	Seat To Floor Height	Safe Working Load	Unit Weight
BTS105300	Aspire Bath Transfer Bench - Bariatric	584 mm	670 mm	390 - 520 mm	227 kg	7.2 kg