

Malvestio Delta 4 Ward Bed - Delta Drive & Weigh Scale



The DELTA4 Ward Bed is the new generation of beds developed according to the needs of healthcare professionals, with an idea in mind: if the bed allows the operators to work at their best, the patient's well-being follows naturally. The DELTA4 bed becomes an active platform that simplifies and maximises the performances of healthcare professionals.

- **INCENTRO System** - the combined back and pelvic translation nullifies horizontal forces, prevents slipping and keeps the body at a 30° angle, optimal for ventilation.
- **Safe Exit Lights** - 2 x LED lights signal whether the bed is at a safe exit height, promoting safe egress for the patient (orange is unsafe and green is safe).
- **Safe Exit control** moves the platform to the minimum height and lifts the backrest up to facilitate the patient to get up from the bed.
- **Auto Regression** - backrest automatically shifts backwards thus increasing the space on the pelvic area by creating an extra space of 13cm.
- **Biocote** - metal components are epoxy powder coated and finished with an antibacterial additive for improved infection control.
- **DeltaDrive** - electric centre 5th wheel with integrated controls in headboard provides increased manoeuvrability for hospital staff (available on BEB050150 only).
- **Weigh Scale** - integrated system allows for precision weight measurements that can be set to automatically weigh a patient every 24 hours (available on BEB050150 only).
- **4 sector safety sides** with integrated electric controls which allows full autonomy of movement to the patient.
- **CPR control** on both sides allows the bed platform to be brought to the horizontal position immediately accessible.
- **Mattress platform** can be easily removed, all sections are designed for easy

✓ Contract approved

Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click [here](#) for more information

cleaning.

- **Automatic self-brake** once connected to the power net a sound alarm alerts the nursing staff about the unbraked bed and after 20seconds the brake is automatically activated.
- **X-ray tray** together with HPL backrest allows extremely easy x-ray operations to be carried out without moving the patient.

Warranty Information

- 10 Years Frame
- 2 Years Electrics
- 1 Year handset

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
BEB050150	Delta 4 Ward Bed - Delta Drive & Weigh Scale	2200 mm	990 mm	390 - 790 mm	300 kg	125 kg