

Jobst UlcerCare Starter Kit



- Enables effective management of chronic lower leg wounds
- Provides graduated compression of 40 mmHg with a combination of a mild compression liner and an outer stocking
- Outer stocking comes with a zipper on the left or right side for comfortable closure on the opposite side of the wound
- Functional with normal footwear and clothing
- UlcerCare Starter Kits include one outer stocking and two liners

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information



Code	Display Name	Calf/Arm Circumference	Ankle Circumference
PTS640990	Jobst UlcerCare - Starter Kit (Stocking + 2 Liners) - 40mmHg - Right Zip - Sml	29-35cm	18-21cm
PTS641000	Jobst UlcerCare - Starter Kit (Stocking + 2 Liners) - 40mmHg - Right Zip - Med	34-41cm	21-24cm
PTS641010	Jobst UlcerCare - Starter Kit (Stocking + 2 Liners) - 40mmHg - Right Zip - Lge	40-47cm	24-27cm
PTS641020	Jobst UlcerCare - Starter Kit (Stocking + 2 Liners) - 40mmHg - Right Zip - XL	46-54cm	27-30cm
PTS641021	Jobst UlcerCare - Starter Kit (Stocking + 2 Liners) - 40mmHg - Right Zip - XXL	48-56cm	30-33cm
PTS641951	Jobst UlcerCare - Starter Kit (Stocking + 2 Liners) - 40mmHg - Left Zip - Sml	29-35cm	18-21cm
PTS641952	Jobst UlcerCare - Starter Kit (Stocking + 2 Liners) - 40mmHg - Left Zip - Med	34-41cm	21-24cm

Code	Display Name	Calf/Arm Circumference	Ankle Circumference
PTS641953	Jobst UlcerCare - Starter Kit (Stocking + 2 Liners) - 40mmHg - Left Zip - Lge	40-47cm	24-27cm
PTS641954	Jobst UlcerCare - Starter Kit (Stocking + 2 Liners) - 40mmHg - Left Zip - XL	46-54cm	27-30cm
PTS641955	Jobst UlcerCare - Starter Kit (Stocking + 2 Liners) - 40mmHg - Left Zip - XXL	48-56cm	30-33cm