

JAY - FOCUS POINT Backrest



JAY[®]

- Designed for the client who has moderate to high postural needs who needs an easy-to-use back support that provides both posterior pelvic and lateral positioning that will change over time
- The Jay Focus Point is a 3 piece backrest featuring a reinforced contoured aluminium shell, multi-adjustable lateral wings and a pelvic wedge for additional posterior pelvic support
- Featuring back height adjustment of 4"; forward tilt (20°); rear tilt (45°); and multiple optional positioning supports, the JAY Focus Point provides the flexibility to adapt to the changing needs of the client
- The shape of the JAY Focus Point minimises upper limb interference
- All models include: Universal mounting hardware and 3 piece cover

Warranty Information

- 2 Years Backrest
- 6 Months Cover

 Contract approved

Support @ Home

Function: Accessories for wheelchair seating

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_181003192_0105_1_2

Wheelchair Manual – Composite:

Wheelchair Powered –

Composite:05_417_0105_1_2

05_419_0105_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Safe Working Load	Seat Width	Backrest Height	Backrest Width
PTB532686	JAY Focus Point Tall - 14 - 15in	113 kg	14 - 15 in	17 in	14 - 15 in

Code	Display Name	Safe Working Load	Seat Width	Backrest Height	Backrest Width
PTB532690	JAY Focus Point - 16 - 17in	113 kg	16 - 17 in	15 in	16 - 17 in
PTB532710	JAY Focus Point - 18 - 19in	113 kg	18 - 19 in	16 in	18 - 19 in
PTB532724	JAY Focus Point - 20 - 21in	113 kg	20 - 21 in	17 in	20 - 21 in
PTB532725	JAY Focus Point Tall - 20 - 21in	113 kg	20 - 21 in	20 in	20 - 21 in