

## Handi Lift Walk Belt - Miracle Grip



- A rectangular belt made from a strong lightweight nylon material
- Four padded webbing handles and optional slip resistant lining
- Two carers, each using both hands, can assist with lifting and / or walking
- The belt has a double hook and loop tape fastening for holding the belt firmly in place
- Available in four sizes
- Inside is non-slip thanks to Miracle Grip

### Warranty Information

- 1 Year

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Changing body position

**Category:** Prescribed

Click [here](#) for more information

### NDIS

**Item Code:**

05\_123115159\_0105\_1\_2

**Minor Transfer –**

**Composite:**05\_430\_0103\_1\_2

Click [here](#) for more information

| Code      | Display Name                                 | Total Length   |
|-----------|----------------------------------------------|----------------|
| TSE673000 | Handi Lift Walk Belt - Miracle Grip - Small  | 690 -930 mm    |
| TSE673120 | Handi Lift Walk Belt - Miracle Grip - Medium | 1100 mm        |
| TSE672990 | Handi Lift Walk Belt - Miracle Grip - Large  | 900 - 1170 mm  |
| TSE672995 | Handi Lift Walk Belt - Miracle Grip - Extra  | 1080 - 1330 mm |

| Code  | Display Name | Total Length |
|-------|--------------|--------------|
| Large |              |              |