

## Handi Lift Walk Belt - Miracle Grip



- A rectangular belt made from a strong lightweight nylon material
- Four padded webbing handles and optional slip resistant lining
- Two carers, each using both hands, can assist with lifting and / or walking
- The belt has a double hook and loop tape fastening for holding the belt firmly in place
- Available in four sizes
- Inside is non-slip thanks to Miracle Grip

### Warranty Information

- 1 Year

✓ Contracted DVA product

✓ Contract approved

**ndis**

**Item Code:**

05\_123115159\_0105\_1\_2

**Minor Transfer –**

**Composite:**05\_430\_0103\_1\_2

Click [here](#) for more information

Code	Display Name	Total Length
TSE673000	Handi Lift Walk Belt - Miracle Grip - Small	690 -930 mm
TSE673120	Handi Lift Walk Belt - Miracle Grip - Medium	1100 mm
TSE672990	Handi Lift Walk Belt - Miracle Grip - Large	900 - 1170 mm
TSE672995	Handi Lift Walk Belt - Miracle Grip - Extra Large	1080 - 1330 mm