

Jobst UltraSheer Compression Stockings - Women



- JOBST UltraSheer compression stockings combine gradient compression therapy with superb style and elegance
- Their fine, breathable satin-touch yarn provides superb sheerness
- The stylish shimmering appearance supports a glamorous look while helping to manage vascular conditions
- Transparent - Sheerness for beautiful legs
- Breathable - Open loop knitting structure promotes high degree of air permeability
- Shimmering - For a glossy look

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information

Warranty Information

- 1 Year

Code	Display Name	Colour	Calf/Arm Circumference	Waist/Hip Circumference	Knee Circumference	Ankle Circumference
PTS640140	Jobst UltraSheer - Knee - 15-20mmHg - C/Toe - Sml - Blk	Black	28-38cm	-	-	18-21cm
PTS640150	Jobst UltraSheer - Knee - 15-20mmHg - C/Toe - Med - Blk	Black	30-42cm	-	-	21-25cm
PTS640160	Jobst UltraSheer - Knee - 15-20mmHg - C/Toe - Lge - Blk	Black	32-46cm	-	-	25-29cm
PTS640170	Jobst UltraSheer - Knee - 15-20mmHg - C/Toe - XL - Blk	Black	34-50cm	-	-	29-33cm
PTS640610	Jobst UltraSheer - Knee - 20-30mmHg - C/Toe - Sml - Blk	Black	28-38cm	-	-	18-21cm
PTS640620	Jobst UltraSheer - Knee - 20-30mmHg - C/Toe - Med - Blk	Black	30-42cm	-	-	21-25cm

Code	Display Name	Colour	Calf/Arm Circumference	Waist/Hip Circumference	Knee Circumference	Ankle Circumference
PTS640630	Jobst UltraSheer - Knee - 20-30mmHg - C/Toe - Lge - Blk	Black	32-46cm	-	-	25-29cm
PTS640640	Jobst UltraSheer - Knee - 20-30mmHg - C/Toe - XL - Blk	Black	34-50cm	-	-	29-33cm
PTS640650	Jobst UltraSheer - Knee - 20-30mmHg - C/Toe - Sml - Natural	Natural	28-38cm	-	-	18-21cm
PTS640660	Jobst UltraSheer - Knee - 20-30mmHg - C/Toe - Med - Natural	Natural	30-42cm	-	-	21-25cm
PTS640670	Jobst UltraSheer - Knee - 20-30mmHg - C/Toe - Lge - Natural	Natural	32-46cm	-	-	25-29cm
PTS640680	Jobst UltraSheer - Knee - 20-30mmHg - C/Toe - XL - Natural	Natural	34-50cm	-	-	29-33cm
PTS640870	Jobst UltraSheer - Knee - 30-40mmHg - C/Toe - Sml - Blk	Black	28-38cm	-	-	18-21cm
PTS640880	Jobst UltraSheer - Knee - 30-40mmHg - C/Toe - Med - Blk	Black	30-42cm	-	-	21-25cm
PTS640890	Jobst UltraSheer - Knee - 30-40mmHg - C/Toe - Lge - Blk	Black	32-46cm	-	-	25-29cm
PTS640900	Jobst UltraSheer - Knee - 30-40mmHg - C/Toe - XL - Blk	Black	34-50cm	-	-	29-33cm
PTS640910	Jobst UltraSheer - Knee - 30-40mmHg - C/Toe - Sml - Natural	Natural	28-38cm	-	-	18-21cm
PTS640920	Jobst UltraSheer - Knee - 30-40mmHg - C/Toe - Med - Natural	Natural	30-42cm	-	-	21-25cm
PTS640930	Jobst UltraSheer - Knee - 30-40mmHg - C/Toe - Lge -	Natural	32-46cm	-	-	25-29cm

Code	Display Name	Colour	Calf/Arm Circumference	Waist/Hip Circumference	Knee Circumference	Ankle Circumference
	Natural					
PTS640940	Jobst UltraSheer - Knee - 30-40mmHg - C/Toe - XL - Natural	Natural	34-50cm	-	-	29-33cm
PTS640300	Jobst UltraSheer - Thigh - 15-20mmHg - C/Toe - Med - Blk	Black	30-42cm	-	46-70cm	21-25cm
PTS640310	Jobst UltraSheer - Thigh - 15-20mmHg - C/Toe - Lge - Blk	Black	32-46cm	-	54-78cm	25-29cm
PTS640315	Jobst UltraSheer - Thigh - 15-20mmHg - C/Toe - XL - Blk	Black	34-50cm	-	60-81cm	29-33cm
PTS640330	Jobst UltraSheer - Thigh - 15-20mmHg - C/Toe - Sml - Natural	Black	28-38cm	-	40-62cm	18-21cm
PTS640335	Jobst UltraSheer - Thigh - 15-20mmHg - C/Toe - Med - Natural	Natural	30-42cm	-	46-70cm	21-25cm
PTS640350	Jobst UltraSheer - Thigh - 15-20mmHg - C/Toe - Lge - Natural	Natural	32-46cm	-	54-78cm	25-29cm
PTS640360	Jobst UltraSheer - Thigh - 15-20mmHg - C/Toe - XL - Natural	Natural	34-50cm	-	60-81cm	29-33cm
PTS640710	Jobst UltraSheer - Thigh - 20-30mmHg - C/Toe - Sml - Blk	Black	28-38cm	-	40-62cm	18-21cm
PTS640720	Jobst UltraSheer - Thigh - 20-30mmHg - C/Toe - Med - Blk	Black	30-42cm	-	46-70cm	21-25cm
PTS640730	Jobst UltraSheer - Thigh - 20-30mmHg - C/Toe - Lge - Blk	Black	32-46cm	-	54-78cm	25-29cm
PTS640735	Jobst UltraSheer - Thigh - 20-30mmHg - C/Toe - XL - Blk	Black	34-50cm	-	60-81cm	29-33cm
PTS640750	Jobst UltraSheer - Thigh - 20-30mmHg - C/Toe - Med - Natural	Natural	30-42cm	-	46-70cm	21-25cm

Code	Display Name	Colour	Calf/Arm Circumference	Waist/Hip Circumference	Knee Circumference	Ankle Circumference
PTS640760	Jobst UltraSheer - Thigh - 20-30mmHg - C/Toe - Lge - Natural	Natural	32-46cm	-	54-78cm	25-29cm
PTS640761	Jobst UltraSheer - Thigh - 20-30mmHg - C/Toe - XL - Natural	Natural	34-50cm	-	60-81cm	29-33cm
PTS640945	Jobst UltraSheer - Thigh - 30-40mmHg - C/Toe - Sml - Blk	Black	28-38cm	-	40-62cm	18-21cm
PTS640946	Jobst UltraSheer - Thigh - 30-40mmHg - C/Toe - Med - Blk	Black	30-42cm	-	46-70cm	21-25cm
PTS640947	Jobst UltraSheer - Thigh - 30-40mmHg - C/Toe - Lge - Blk	Black	32-46cm	-	54-78cm	25-29cm
PTS640948	Jobst UltraSheer - Thigh - 30-40mmHg - C/Toe - XL - Blk	Black	34-50cm	-	60-81cm	29-33cm
PTS640949	Jobst UltraSheer - Thigh - 30-40mmHg - C/Toe - Sml - Natural	Natural	28-38cm	-	40-62cm	18-21cm
PTS640950	Jobst UltraSheer - Thigh - 30-40mmHg - C/Toe - Med - Natural	Natural	30-42cm	-	46-70cm	21-25cm
PTS640955	Jobst UltraSheer - Thigh - 30-40mmHg - C/Toe - Lge - Natural	Natural	32-46cm	-	54-78cm	25-29cm
PTS640960	Jobst UltraSheer - Thigh - 30-40mmHg - C/Toe - XL - Natural	Natural	34-50cm	-	60-81cm	29-33cm
PTS640380	Jobst UltraSheer - Waist - 15-20mmHg - C/Toe - Sml - Blk	Black	28-38cm	71-117cm	40-62cm	18-21cm
PTS640390	Jobst UltraSheer - Waist - 15-20mmHg - C/Toe - Med - Blk	Black	30-42cm	76-127cm	46-70cm	21-25cm
PTS640400	Jobst UltraSheer - Waist - 15-20mmHg - C/Toe - Lge - Blk	Black	32-46cm	81-137cm	54-78cm	25-29cm
PTS640410	Jobst UltraSheer - Waist - 15-20mmHg - C/Toe - XL - Blk	Black	34-50cm	102-166cm	60-81cm	29-33cm

Code	Display Name	Colour	Calf/Arm Circumference	Waist/Hip Circumference	Knee Circumference	Ankle Circumference
PTS640450	Jobst UltraSheer - Waist - 15-20mmHg - C/Toe - Sml - Natural	Natural	28-38cm	71-117cm	40-62cm	18-21cm
PTS640460	Jobst UltraSheer - Waist - 15-20mmHg - C/Toe - Med - Natural	Natural	30-42cm	76-127cm	46-70cm	21-25cm
PTS640470	Jobst UltraSheer - Waist - 15-20mmHg - C/Toe - Lge - Natural	Natural	32-46cm	81-137cm	54-78cm	25-29cm
PTS640480	Jobst UltraSheer - Waist - 15-20mmHg - C/Toe - XL - Natural	Natural	34-50cm	102-166cm	60-81cm	29-33cm
PTS640779	Jobst UltraSheer - Waist - 20-30mmHg - C/Toe - Sml - Blk	Black	28-38cm	71-117cm	40-62cm	18-21cm
PTS640780	Jobst UltraSheer - Waist - 20-30mmHg - C/Toe - Med - Blk	Black	30-42cm	76-127cm	46-70cm	21-25cm
PTS640781	Jobst UltraSheer - Waist - 20-30mmHg - C/Toe - Lge - Blk	Black	32-46cm	81-137cm	54-78cm	25-29cm
PTS640782	Jobst UltraSheer - Waist - 20-30mmHg - C/Toe - XL - Blk	Black	34-50cm	102-166cm	60-81cm	29-33cm
PTS640789	Jobst UltraSheer - Waist - 20-30mmHg - C/Toe - Sml - Natural	Natural	28-38cm	71-117cm	40-62cm	18-21cm
PTS640790	Jobst UltraSheer - Waist - 20-30mmHg - C/Toe - Med - Natural	Natural	30-42cm	76-127cm	46-70cm	21-25cm
PTS640800	Jobst UltraSheer - Waist - 20-30mmHg - C/Toe - Lge - Natural	Natural	32-46cm	81-137cm	54-78cm	25-29cm
PTS640810	Jobst UltraSheer - Waist - 20-30mmHg - C/Toe - XL - Natural	Natural	34-50cm	102-166cm	60-81cm	29-33cm
PTS640969	Jobst UltraSheer - Waist - 30-40mmHg - C/Toe - Sml -	Natural	28-38cm	71-117cm	40-62cm	18-21cm

Code	Display Name	Colour	Calf/Arm Circumference	Waist/Hip Circumference	Knee Circumference	Ankle Circumference
	Natural					
PTS640970	Jobst UltraSheer - Waist - 30-40mmHg - C/Toe - Med - Natural	Natural	30-42cm	76-127cm	46-70cm	21-25cm
PTS640980	Jobst UltraSheer - Waist - 30-40mmHg - C/Toe - Lge - Natural	Natural	32-46cm	81-137cm	54-78cm	25-29cm
PTS640981	Jobst UltraSheer - Waist - 30-40mmHg - C/Toe - XL - Natural	Natural	34-50cm	102-166cm	60-81cm	29-33cm
PTS640983	Jobst UltraSheer - Waist - 30-40mmHg - C/Toe - Sml - Blk	Black	28-38cm	71-117cm	40-62cm	18-21cm
PTS640984	Jobst UltraSheer - Waist - 30-40mmHg - C/Toe - Med - Blk	Black	30-42cm	76-127cm	46-70cm	21-25cm
PTS640986	Jobst UltraSheer - Waist - 30-40mmHg - C/Toe - XL - Blk	Black	34-50cm	102-166cm	60-81cm	29-33cm
PTS647015	JOBST UltraSheer - Knee - 15-20mmHg - C/Toe - Lge - Natural	Natural	-	-	-	-
PTS647020	JOBST UltraSheer - Knee - 15-20mmHg - C/Toe - Med - Natural	Natural	-	-	-	-
PTS647025	JOBST UltraSheer - Knee - 15-20mmHg - C/Toe - Sml - Natural	Natural	-	-	-	-
PTS647030	JOBST UltraSheer - Knee - 15-20mmHg - C/Toe - XL - Natural	Natural	-	-	-	-
PTS640220	JOBST UltraSheer - Knee - 15-20mmHg - C/Toe - Sml - Natural	Natural	-	-	-	-
PTS640230	JOBST UltraSheer - Knee - 15-20mmHg - C/Toe - Med -	Natural	-	-	-	-

Code	Display Name	Colour	Calf/Arm Circumference	Waist/Hip Circumference	Knee Circumference	Ankle Circumference
	Natural					
PTS640240	JOBST UltraSheer - Knee - 15-20mmHg - C/Toe - Lge - Natural	Natural -	-	-	-	-
PTS640250	JOBST UltraSheer - Knee - 15-20mmHg - C/Toe - XL - Natural	Natural -	-	-	-	-
PTS644910	Jobst Ultrasheer - Thigh W/ Grip Top - 15-20Mmhg - C/Toe - SML - Blk	Mexico -	-	-	-	-